



SPECIFICATION

Benzamidine HCl

Product Code:	FB06195
CAS Number:	1670-14-0
Chemical Formula:	$C_7H_8N_2 \cdot HCl$
Molecular Weight:	156.61

TECHNICAL SPECIFICATION

Appearance:	White crystals
Assay (HPLC) (area %):	min 99%
Assay (Titration - chloride):	min 98%
Total impurities (HPLC):	max 2%
Single impurity (HPLC):	max 0.1%
Water content (Karl Fischer):	max 1%