



SPECIFICATION

UDP-N-acetyl-D-galactosamine

Product Code:	MU04515
Synonyms:	2-Acetamido-2-deoxy-a-D-galactopyranosyl uridine diphosphate UDP-GalNAc Uridine 5'-diphospho-N-acetylgalactosamine disodium salt
CAS Number:	108320-87-2
Chemical Formula:	$C_{17}H_{25}N_3Na_2O_{17}P_2$
Molecular Weight:	651.32

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 95%
Water content (Karl Fischer):	max 7%