



SPECIFICATION

D-Lyxono-1,4-lactone

| | |
|-------------------|-------------|
| Product Code: | ML05267 |
| CAS Number: | 15384-34-6 |
| Chemical Formula: | $C_5H_8O_5$ |
| Molecular Weight: | 148.11 |

TECHNICAL SPECIFICATION

| | | |
|---|---------------------------------------|--|
| Appearance: | White to off-white crystalline powder | |
| Purity (HPLC): | min 97% | |
| Specific optical rotation $[\alpha]_D$ | +80 to +90° | ²⁰ (c=1, H ₂ O, 24 hrs): |