



SPECIFICATION

TBTU

Product Code:	FT07096
Synonyms:	N,N,N',N'-Tetramethyl-O-(benzotriazol-1-yl)uronium tetrafluoroborate
CAS Number:	125700-67-6
Chemical Formula:	$C_{11}H_{16}N_5O.BF_4$
Molecular Weight:	321.08

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 99%
Loss on drying:	max 0.5%
Melting point:	min 200°C (dec.)