



SPECIFICATION

Lewis B tetrasaccharide

| | |
|-------------------|---------------------------------------|
| Product Code: | OL02434 |
| Synonyms: | Fuc-a-1,2-Gal-b-1,3-(Fuc-a-1,4)GlcNAc |
| CAS Number: | 80081-06-7 |
| Chemical Formula: | $C_{26}H_{45}NO_{19}$ |
| Molecular Weight: | 675.63 |

TECHNICAL SPECIFICATION

| | |
|------------------------|-----------------------|
| Appearance: | White powder |
| Purity (TLC): | min 95% |
| Identity (1H NMR): | Conforms to structure |