



SPECIFICATION

Telmisartan acyl-b-D-glucuronide

Product Code:	MT07380
Synonyms:	1-[4'-[[1',4-Dimethyl-2-propyl[2'-propyl[2',6-b1-1H-benzimidazol]-1-yl)methyl]][1,1'-biphenyl]-2-carboxylate]-b-D-glucopyranuronic acid
CAS Number:	250780-40-6
Chemical Formula:	$C_{39}H_{38}N_4O_8$
Molecular Weight:	690.74

TECHNICAL SPECIFICATION

Appearance:	Off-white crystalline powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure