



SPECIFICATION

Lewis A

Product Code:	OL04541
Synonyms:	2-Acetamido-2-deoxy-4-O-(α -L-fucopyranosyl)-3-O-(β -D-glucopyranosyl)-D-glucopyranose Gal- β -1,3(Fuc- α -1,4)GlcNAc
CAS Number:	56570-03-7
Chemical Formula:	$C_{20}H_{35}NO_{15}$
Molecular Weight:	529.49

TECHNICAL SPECIFICATION

Appearance:	White lyophilised solid
Purity (TLC):	min 95%
Identity (NMR) (D ₂ O):	Conforms to structure