



SPECIFICATION

Luteolin-7-O-D-glucopyranoside

Product Code:	ML08112
Synonyms:	2-(3,4-Dihydroxyphenyl)-5-hydroxy-4-oxo-4H-chromen-7-yl-D-glucopyranoside Cynaroside Luteoloside
CAS Number:	5373-11-5
Chemical Formula:	$C_{21}H_{20}O_{11}$
Molecular Weight:	448.38

TECHNICAL SPECIFICATION

Appearance:	Yellow powder
Assay (HPLC):	min 98%