



SPECIFICATION

5-Bromo-4-chloro-3-indolyl b-D-glucuronide sodium salt

Product Code:	EB03265
Synonyms:	X-b-D-Glucuronide X-b-D-GLcA-Na
CAS Number:	370100-64-4
Chemical Formula:	$C_{14}H_{12}BrClNaO_7$
Molecular Weight:	444.59

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 98%
Solubility (1% in water):	Clear and practically colourless solution
Water content (Karl Fischer):	max 10%