



SPECIFICATION

6-Bromo-2-naphthyl b-D-galactopyranoside

| | |
|-------------------|---------------------|
| Product Code: | EB04969 |
| CAS Number: | 15572-30-2 |
| Chemical Formula: | $C_{16}H_{17}BrO_6$ |
| Molecular Weight: | 385.21 |

TECHNICAL SPECIFICATION

| | |
|-----------------|---|
| Appearance: | White or slightly yellow crystalline powder |
| Purity (HPLC): | min 98% |
| Loss on Drying: | max 1.0% |
| Identity (IR): | Conforms to structure |
| Heavy Metals: | max 20ppm |