



SPECIFICATION

2-Acetamido-3,4,6-tri-O-acetyl-2-deoxy-b-D-glucopyranosyl azide

Product Code:	MA00690
CAS Number:	6205-69-2
Chemical Formula:	$C_{14}H_{20}N_4O_8$
Molecular Weight:	372.33

TECHNICAL SPECIFICATION

Appearance:	White crystalline powder
Purity (HPLC):	min 98%
Loss on Drying:	max 1.0%
Melting Point:	163 - 166°C
Identity (IR):	Conforms to structure
Solubility:	Easily soluble in EtOAc, insoluble in H ₂ O