



SPECIFICATION

N4-Acetyl-2'-deoxy-5'-O-DMT-2'-fluorocytidine

Product Code:	NA06052
Synonyms:	2'-Fluoro-5'-O-DMT-N4-acetyl-cytidine N4-Acetyl-5'-O-DMT-2'-fluoro-2'-deoxycytidine
CAS Number:	159414-98-9
Chemical Formula:	$C_{32}H_{32}FN_3O_7$
Molecular Weight:	589.61

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 98%
Residual solvents:	max 5%