



SPECIFICATION

Lewis X trisaccharide

Product Code:	OL06490
Synonyms:	Gal-b-1,4(Fuc-a1,3)GlcNAc 2-Acetamido-2-deoxy-3-O-(a-L-fucopyranosyl)-4-O-(b-D-galactopyranosyl)-D-glucopyranoside
CAS Number:	71208-06-5
Chemical Formula:	$C_{20}H_{35}NO_{15}$
Molecular Weight:	529.49

TECHNICAL SPECIFICATION

Appearance:	White lyophilised solid
Purity (TLC):	min 95%
Identity (NMR):	Conforms to structure