



SPECIFICATION

Guanosine 5-O-(2-thiodiphosphate), trilithium Salt

Product Code:	NG09615
CAS Number:	97952-36-8
Chemical Formula:	$C_{10}H_{12}N_5O_{10}P_2S_2 \cdot 3Li$
Molecular Weight:	477.1

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 85%
Impurities (GMP):	max 10%
Lithium content:	max 5%
Water content (Karl Fischer):	max 12%