



SPECIFICATION

S-Adenosyl-L-methionine 1,4-butanedisulfonate sodium salt

Product Code:	NA08265
CAS Number:	101020-79-5
Chemical Formula:	$C_{15}H_{24}N_6O_5S.C_4H_8Na_2O_8S_2$
Molecular Weight:	694.67

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Identity (IR):	Conforms to structure
Heavy metals:	max 10ppm
Butanedisulfonic acid (Titration):	46 -49%
Methionine (HPLC):	50 - 53%