



## SPECIFICATION

### N-(5-Carboxypentyl)-deoxynojirmycin

Product Code:	MC04528
Synonyms:	N-(5-Carboxypentyl)-1,5-dideoxy-1,5-imino-D-glucitol
CAS Number:	79206-70-5, 79206-51-2
Chemical Formula:	$C_{12}H_{23}NO_6$
Molecular Weight:	277.31

## TECHNICAL SPECIFICATION

Appearance:	White to beige powder
Purity (TLC):	min 98%
Identity ( $^1H$ NMR):	Conforms to structure