



SPECIFICATION

3,4,6-Tri-O-acetyl-D-galactal

Product Code:	MT05853
CAS Number:	4098-06-0
Chemical Formula:	$C_{12}H_{16}O_7$
Molecular Weight:	272.25

TECHNICAL SPECIFICATION

Appearance:	Light yellow solidifying oil
Purity (GC):	min 98%
Specific optical rotation $[\alpha]_D^{20}$	-17 to -22° ²⁰ (c=3, CHCl ₃):
Identity (¹ H NMR):	Conforms to structure
Residual solvents:	max 5%