



SPECIFICATION

5-Bromo-4-chloro-3-indolyl b-D-glucuronide sodium salt trihydrate

Product Code:	EB03171
Synonyms:	X-Gluc sodium salt
CAS Number:	129541-41-9
Chemical Formula:	$C_{14}H_{12}BrClNaO_7 \cdot 3H_2O$
Molecular Weight:	498.64

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 98%
Solubility 1% w/v:	Clear and almost colourless