



SPECIFICATION

Voglibose

Product Code:	MV06910
Synonyms:	3,4-Dideoxy-[2-hydroxy-1-(hydroxyethyl)ethyl] amino-1-C-(hydroxynethyl)-1,2,3,4-cyclohexaneterol 3,4-Dideoxy-4-[[2-hydroxy-1-(hydroxymethyl)ethyl]amino]-2-C-(hydroxymethyl)-D-epiinositol
CAS Number:	83480-29-9
Chemical Formula:	$C_{10}H_{21}NO_6$
Molecular Weight:	251.28

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Assay:	min 98%
Residue on ignition:	max 0.1%
Water content (Karl Fischer):	max 1%