



SPECIFICATION

S-(4-Nitrobenzyl)-6-thioinosine

Product Code:	NN09410
Synonyms:	6-[(4-Nitrobenzyl)thio]-9-b-D-ribofuranosylpurine
CAS Number:	38048-32-7
Chemical Formula:	$C_{17}H_{17}N_5O_6S$
Molecular Weight:	419.41

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity:	min 98%
Identity (1H NMR):	Conforms to structure