



## SPECIFICATION

### 6-Bromo-2-naphthyl b-D-glucopyranoside

Product Code:	EB03284
CAS Number:	15548-61-5
Chemical Formula:	$C_{16}H_{17}BrO_6$
Molecular Weight:	385.21

#### TECHNICAL SPECIFICATION

Appearance:	White or slightly yellow crystalline powder
Purity (HPLC):	min 98%
Loss on Drying:	max 1.0%
Identity (IR):	Conforms to structure
Heavy Metals:	max 20ppm