



SPECIFICATION

Sodium-D-glucuronate

Product Code:	MS04198
Synonyms:	D-Glucuronic acid sodium salt
CAS Number:	14984-34-0, 7182-77-6
Chemical Formula:	$C_6H_9NaO_7$
Molecular Weight:	216.12

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (Titration):	min 98%
Specific optical rotation $[\alpha]_D$	+21 to +23° ²⁰ (c=2, H ₂ O):
Solubility (50 mg/ml in H ₂ O):	Clear colourless to yellow solution
Water content (Karl Fischer):	max 8.5%