



SPECIFICATION

5-Bromo-4-chloro-3-indolyl b-D-glucopyranoside

Product Code:	EB03169
Synonyms:	X-Glu X-b-D-Glucoside
CAS Number:	15548-60-4
Chemical Formula:	$C_{14}H_{15}BrClNO_6$
Molecular Weight:	408.63

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder	
Purity (HPLC):	min 98%	
Specific Optical Rotation $[\alpha]_D$	$-81^\circ \pm 2^\circ$	²⁰ (c=1, DMF):
Identity (IR):	Conforms to structure	
Water content (Karl Fischer):	max 1%	