



## SPECIFICATION

### 3-Indolyl b-D-galactopyranoside

Product Code:	EI05985
Synonyms:	3-Indoxyl b-D-galactopyranoside
CAS Number:	126787-65-3
Chemical Formula:	$C_{14}H_{17}NO_6$
Molecular Weight:	295.29

## TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 98%
Specific optical rotation $[\alpha]_D$	-56 to -52° <sup>20</sup> (c=0.5, MeOH):
Identity (IR):	Conforms to structure