



## SPECIFICATION

### 3-O-(2-Acetamido-2-deoxy- $\alpha$ -galactopyranosyl)-D-galactopyranoside

Product Code:	OA08085
Synonyms:	$\alpha$ -D-N-Acetylgalactosaminyl 1,3 galactose
Chemical Formula:	$C_{14}H_{25}NO_{11}$
Molecular Weight:	383.33

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (TLC):	min 95%
Identity ( $^1H$ NMR):	Conforms to structure