



SPECIFICATION

Sedoheptulose anhydride

Product Code:	MS08092
Synonyms:	2,7-Anhydro-b-D-altro-heptulopyranose
CAS Number:	469-90-9
Chemical Formula:	$C_7H_{12}O_6$
Molecular Weight:	192.17

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 90%
Specific optical rotation $[\alpha]_D$	-129 to -133° ²⁰ (c=2, H ₂ O):