



## SPECIFICATION

### 1,2,3,6-Tetrahydropyridine

Product Code:	FT06292
CAS Number:	694-05-3
Chemical Formula:	C <sub>5</sub> H <sub>9</sub> N
Molecular Weight:	83.13

## TECHNICAL SPECIFICATION

Appearance:	Clear to slightly yellow liquid
Purity (GC):	min 97%
Identity ( <sup>1</sup> H NMR) (CDCl <sub>3</sub> ):	Conforms to structure
Refractive index n <sub>D</sub> <sup>20</sup> :	1.4800 ± 0.003
Water content (Karl Fischer):	max 1%