



SPECIFICATION

Lacto-N-fucopentaose III

Product Code:	OL04212
Synonyms:	Lewis X pentasaccharide LNFIII Gal-b-1,4-(Fuc-a-1,3)-GlcNAc-b-1,3-D-Gal-b-1,4-Glc
CAS Number:	25541-09-7
Chemical Formula:	$C_{32}H_{55}NO_{25} \cdot H_2O$
Molecular Weight:	871.79

TECHNICAL SPECIFICATION

Appearance:	White freeze-dried powder
Purity (HPLC):	min 95%