



## SPECIFICATION

### 2-Acetamido-1,3,4,6-tetra-O-acetyl-2-deoxy- $\alpha$ -D-galactopyranose

Product Code:	MA08854
CAS Number:	10385-50-9
Chemical Formula:	$C_{16}H_{23}NO_{10}$
Molecular Weight:	389.36

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Assay:	min 98%
Purity (TLC):	One spot
Specific optical rotation [ $\alpha$ ] <sub>D</sub>	+99 to +105° <sup>20</sup> (c=1, CHCl <sub>3</sub> ):
Identity ( <sup>1</sup> H NMR):	Conforms to structure