



SPECIFICATION

Vindoline

Product Code:	FV10830
Synonyms:	(2b,3b,4b,5a,12b,19a)-4-(Acetyloxy)-6,7-didehydro-3-hydroxy-16-methoxy-1-methyl-aspidospermidine-3-carboxylic Acid Methyl Ester;(-)-Vindoline NSC 91994
CAS Number:	2182-14-1
Chemical Formula:	$C_{25}H_{32}N_2O_6$
Molecular Weight:	456.53

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 98%
Identity (IR):	Conforms to structure
Loss on drying:	max 1.0%
Residue on ignition:	max 0.2%