



## SPECIFICATION

### Raloxifene 6'-D-glucuronide

Product Code:	MR07278
Synonyms:	2-(4-Hydroxyphenyl)-3-[4-2[2-(1-piperindinyl)ethoxybenzoyl]benzo[b]thien-6-yl]-b-D-glucopyranosiduronic acid
CAS Number:	174264-50-7
Chemical Formula:	C <sub>34</sub> H <sub>35</sub> NO <sub>10</sub> S
Molecular Weight:	649.71

## TECHNICAL SPECIFICATION

Appearance:	Yellow crystalline powder
Purity (TLC):	min 98%
Identity ( <sup>1</sup> H NMR):	Conforms to structure