



SPECIFICATION

Octadecylthioethyl 4-O-(4-O[6-O-a-D-glucopyranosyl-a-D-glucopyranosyl]-a-D-glucopyranosyl)-b-D-glucopyranoside

Product Code:	DO05214
CAS Number:	90215-01-3
Chemical Formula:	$C_{44}H_{82}O_{21}S$
Molecular Weight:	979.18

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 95%
Water content (Karl Fischer):	max 5%