



SPECIFICATION

2-Acetamido-2-deoxy-4-O-(α -L-fucopyranosyl)-D-glucopyranose

Product Code:	OA04447
Synonyms:	Fuc-a-1,4-GlcNAc
CAS Number:	76211-71-7
Chemical Formula:	$C_{14}H_{25}NO_{10}$
Molecular Weight:	367.35

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure