



## SPECIFICATION

### 2-Acetamido-2-deoxy-4-O-(6-sulfo-b-D-galactopyranosyl)-b-D-glucopyranose sodium salt

Product Code:	OA08362
Synonyms:	N-Acetyllactosamine 6'-sulfate sodium salt
Chemical Formula:	$C_{14}H_{24}NNaO_{14}S$
Molecular Weight:	485.39

## TECHNICAL SPECIFICATION

Appearance:	White to off-white solid
Purity (TLC):	min 95%
Identity (NMR) (D <sub>2</sub> O):	Conforms to structure