



SPECIFICATION

5-Bromo-3-indolyl-b-D-galactopyranoside-6-sulfate sodium salt

Product Code:	EB03281
Chemical Formula:	$C_{14}H_{15}BrNNaO_9S$
Molecular Weight:	476.23

TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder	
Purity (HPLC):	min 98%	
Specific optical rotation $[\alpha]_D$	$-24.5^\circ \pm 2$	$^{20} (c=1, H_2O):$
Water content (Karl Fischer):	max 4.8%	