



SPECIFICATION

4-Chloro-3-indolyl b-D-galactopyranoside

Product Code:	EC06387
CAS Number:	135313-63-2
Chemical Formula:	$C_{14}H_{16}ClNO_6$
Molecular Weight:	329.73

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure
Solubility (50mg/ml in MeOH):	Clear, colourless to pale yellow solution
Water content (Karl Fischer):	Reported