



## SPECIFICATION

### 2-Acetamido-1,2-dideoxynojirimycin

Product Code:	MA00005
Synonyms:	2-Acetamido-1,5-imino-1,2,5-trideoxy-D-glucitol AdNJ N-[(3S,4R,5R,6R)-4,5-Dihydroxy-6-(hydroxymethyl)-3-piperidyl]acetamide
CAS Number:	105265-96-1
Chemical Formula:	$C_8H_{16}N_2O_4$
Molecular Weight:	204.22

## TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 98%
Identity ( $^1H$ NMR):	Conforms to structure