



## SPECIFICATION

### Uracil

Product Code:	FU01843
Synonyms:	2,4(1H,3H)-Pyrimidinedione 2,4-Dihydropyrimidine 2,4-Pyrimidinediol
CAS Number:	66-22-8
Chemical Formula:	$C_4H_4N_2O_2$
Molecular Weight:	112.09

## TECHNICAL SPECIFICATION

Appearance:	White to off-white crystalline powder
Purity (HPLC):	min 98%
Loss on Drying:	max 0.5%
Residue on Ignition:	max 0.1%
Heavy Metals:	max 20ppm