



SPECIFICATION

Resorufin b-D-glucuronide sodium salt

Product Code:	ER03328
CAS Number:	125440-91-7
Chemical Formula:	$C_{18}H_{14}NNaO_9$
Molecular Weight:	411.29

TECHNICAL SPECIFICATION

Appearance:	orange to red-brown powder
Purity (HPLC):	min 90%
Identity (IR):	Conforms to structure
Impurity ((Resorufin)) (HPLC):	max 0.5%
Water content (Karl Fischer):	max 5%
Solubility (0.05% in water):	Clear to slightly hazy orange liquid