



SPECIFICATION

Diphenylcyclopropenone

| | |
|-------------------|-----------------|
| Product Code: | FD06273 |
| CAS Number: | 886-38-4 |
| Chemical Formula: | $C_{15}H_{10}O$ |
| Molecular Weight: | 206.24 |

TECHNICAL SPECIFICATION

| | |
|--|-----------------------------|
| Appearance: | White to light brown powder |
| Purity (HPLC): | min 99% |
| Melting Point: | 119 - 122°C |
| Impurities (TLC) (any single impurity): | max 1% |
| Solubility in water: | Insoluble |