



## SPECIFICATION

### 2-Iodo-2',3'-O-isopropylideneadenosine

Product Code:	NI04752
Synonyms:	2-Iodo-D-adenosine 2',3'-acetonide
CAS Number:	141018-25-9
Chemical Formula:	$C_{13}H_{16}IN_5O_4$
Molecular Weight:	433.21

#### TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	min 98%
Identity ( $^1H$ NMR):	Conforms to structure