



## SPECIFICATION

### D-Maltitol

|                   |                      |
|-------------------|----------------------|
| Product Code:     | OM06867              |
| CAS Number:       | 585-88-6             |
| Chemical Formula: | $C_{12}H_{24}O_{11}$ |
| Molecular Weight: | 344.31               |

## TECHNICAL SPECIFICATION

|  |  |
|--|--|
| Appearance:                                    | White to off-white crystalline powder          |
| Purity (HPLC):                                 | min 99%  |
| Specific Optical Rotation<br>$[\alpha]_D^{20}$ | +104 to +109° $^{20}$ (c=1, H <sub>2</sub> O): |
| Water content (Karl Fischer):                  | max 0.5%                                       |