



SPECIFICATION

Maltotetraose

Product Code:	OM02796
Synonyms:	α -D-Glucopyranosyl-1,4-O- α -D-glucopyranosyl-1,4-O- α -D-glucopyranosyl-1,4-D-glucopyranose
CAS Number:	34612-38-9
Chemical Formula:	$C_{24}H_{42}O_{21}$
Molecular Weight:	666.58

TECHNICAL SPECIFICATION

Appearance:	White powder	
Purity:	min 97%	
Specific optical rotation $[\alpha]_D$	$+177 \pm 1^\circ$	²⁰ (c=3, H ₂ O):
Moisture:	max 5%	