



SPECIFICATION

Dodecyl b-D-maltopyranoside

Product Code:	DD06199
Synonyms:	Dodecyl 4-O-(α -D-glucopyranosyl)-b-D-glucopyranoside
CAS Number:	69227-93-6
Chemical Formula:	$C_{24}H_{46}O_{11}$
Molecular Weight:	510.62

TECHNICAL SPECIFICATION

Appearance:	White powder
Purity (HPLC):	min 99%
Specific optical rotation $[\alpha]_D$	+45 to +49° ²⁰ (c=1, H ₂ O):