



SPECIFICATION

2,2,6,6-Tetramethylpiperidine 1-oxyl, free radical free flowing milled solid

Product Code:	FT09614
Synonyms:	TEMPO
CAS Number:	2564-83-2
Chemical Formula:	$C_9H_{18}NO$
Molecular Weight:	156.25

TECHNICAL SPECIFICATION

Appearance:	Orange to red solid
Assay (GC):	Min 98%
Identity (IR):	Conforms to structure
Residual solvents (as methanol):	Max 0.3%