



SPECIFICATION

Epigallocatechin gallate

Product Code:	FE09970
Synonyms:	(-)-cis-2-(3,4,5-Trihydroxyphenyl)-3,4-dihydro-1(2H)-benzopyran-3,5,7-triol 3-gallate (-)-cis-3,3,4,5,5,7-Hexahydroxy-flavane-3-gallate EGCG
CAS Number:	989-51-5
Chemical Formula:	$C_{22}H_{18}O_{10}$
Molecular Weight:	442.37

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (HPLC):	min 98%
Sulfated ash:	max 0.1%
Water content:	max 0.5%