



SPECIFICATION

Eriocitrin

Product Code:	OE08285
Synonyms:	(S)-3',4',5,7-Tetrahydroxyflavanone-7-[6-O-(α -L-rhamnopyranosyl)- β -D-glucopyranoside] Eriodictyol 7-O-rutinoside
CAS Number:	13463-28-0
Chemical Formula:	$C_{27}H_{32}O_{15}$
Molecular Weight:	596.53

TECHNICAL SPECIFICATION

Appearance:	White to yellow powder
Purity (HPLC):	min 98%
Identity Mass spec:	Conforms to structure
Water content (Karl Fischer):	max 3%