



SPECIFICATION

1-Phenylsulfonylpyrrole-3-carboxylic acid

Product Code:	FP07112
Synonyms:	N-Tosyl-3-pyrrolecarboxylic acid
CAS Number:	106058-86-0
Chemical Formula:	$C_{11}H_9NO_4S$
Molecular Weight:	251.26

TECHNICAL SPECIFICATION

Appearance:	Light brown powder
Purity (HPLC):	min 98%
Melting Point:	202 - 206°C
Identity (1H NMR):	Conforms to structure