



SPECIFICATION

4-Nitrophenyl 2-acetamido-3-O-(2-acetamido-2-deoxy-b-D-glucopyranosyl)-2-deoxy-a-D-galactopyranoside

Product Code:	EN05763
Synonyms:	GlcNAc-b-1,3-GalNAc-a-PNP
CAS Number:	125455-64-3
Chemical Formula:	$C_{22}H_{31}N_3O_{13}$
Molecular Weight:	545.49

TECHNICAL SPECIFICATION

Appearance:	White to off-white powder
Purity (TLC):	Min 98%
Identity (1H NMR):	Conforms to structure