



SPECIFICATION

4-Nitrophenyl 2-acetamido-2-deoxy-3-O-(b-D-galactopyranosyl)-b-D-glucopyranoside

Product Code:	EN06439
Synonyms:	Gal-b-1-3-GlcNAc-b-PNP
CAS Number:	57467-13-7
Chemical Formula:	$C_{20}H_{28}N_2O_{13}$
Molecular Weight:	504.44

TECHNICAL SPECIFICATION

Appearance:	White to light yellow crystalline powder
Purity (TLC):	min 98%
Identity (1H NMR):	Conforms to structure